

— ABLYBIO, Help Your Research



SLC9A3R2 Antibody

货号: **AYP6154**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WBELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	apical plasma membrane,extracellular exosome,focal adhesion,nucleus,plasma membrane
纯化	Affinity purification

抗原信息

抗原信息 Synthesized peptide derived from Human SLC9A3R2.

靶点信息

研究背景 This gene encodes a member of the NHERF family of PDZ scaffolding proteins. These proteins mediate many cellular processes by binding to and regulating the membrane expression and protein-protein interactions of membrane receptors and transport proteins. The encoded protein plays a role in intestinal sodium absorption by regulating the activity of the sodium/hydrogen exchanger 3, and may also regulate the cystic fibrosis transmembrane regulator (CFTR) ion channel. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011]

基因ID	9351
基因名	SLC9A3R2
Swiss	Q15599
别名	E3KARP,NHE3RF2,NHERF-2,NHERF2,OCTS2,SIP-1,SIP1,TKA-1,SLC9A3R2,SLC9A3R2 Antibody,NHE3 kinase A regulatory protein E3KARP,SRY-interacting protein 1,Sodium-hydrogen exchanger regulatory factor 2,Solute carrier family 9 isoform A3 regulatory factor 2,Tyrosine kinase activator protein 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn