

— ABLYBIO, Help Your Research



NDUFA9 Antibody

货号: **AYP6155**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	42kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,BT-474,HL-60,HT-29,Mouse kidney,Mouse heart,Mouse brain,Mouse liver,Rat testis
细胞定位	Mitochondrion matrix
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human NDUFA9.
------	--

靶点信息

研究背景	The encoded protein is a subunit of the hydrophobic protein fraction of the NADH:ubiquinone oxidoreductase (complex I), the first enzyme complex in the electron transport chain located in the inner mitochondrial membrane. A pseudogene has been identified on chromosome 12.
基因ID	4704
基因名	NDUFA9
Swiss	Q16795 (https://www.uniprot.org/uniprotkb/Q16795/entry)
别名	NDUFA9,CC6,C1-39k,C139k,NDUFS2L,SDR22E1,NDUFA9 Antibody,Complex I-39kD,NADH-ubiquinone oxidoreductase 39 kDa subunit

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)