

— ABLYBIO, Help Your Research



# NUP107 Antibody

货号: **AYP6169**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	79kDa/103kDa/106kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,HeLa,293T,Jurkat,HepG2,Mouse brain,Mouse kidney
细胞定位	Chromosome,Nucleus,Nucleus membrane,centromere,kinetochore,nuclear pore complex
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human NUP107.
------	--

## 靶点信息

研究背景	This gene encodes a member of the nucleoporin family. The protein is localized to the nuclear rim and is an essential component of the nuclear pore complex (NPC). All molecules entering or leaving the nucleus either diffuse through or are actively transported by the NPC. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized.
基因ID	57122
基因名	NUP107
Swiss	P57740 ( <a href="https://www.uniprot.org/uniprotkb/P57740/entry">https://www.uniprot.org/uniprotkb/P57740/entry</a> )
别名	NUP107,NPHS11,NUP84,NUP107 Antibody,107 kDa nucleoporin,Nucleoporin Nup107

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)