

— ABLYBIO, Help Your Research



# SLC15A1/PEPT1 Antibody

货号: **AYP6173**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	78kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Rat testis
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human SLC15A1/PEPT1.
------	---

## 靶点信息

研究背景	This gene encodes an intestinal hydrogen peptide cotransporter that is a member of the solute carrier family 15. The encoded protein is localized to the brush border membrane of the intestinal epithelium and mediates the uptake of di- and tripeptides from the lumen into the enterocytes. This protein plays an important role in the uptake and digestion of dietary proteins. This protein also facilitates the absorption of numerous peptidomimetic drugs.
基因ID	6564
基因名	SLC15A1
Swiss	P46059 ( <a href="https://www.uniprot.org/uniprotkb/P46059/entry">https://www.uniprot.org/uniprotkb/P46059/entry</a> )
别名	SLC15A1,HPECT1,HPEPT1,PEPT1,SLC15A1/PEPT1 Antibody,Intestinal H(+)/peptide cotransporter,Oligopeptide transporter,small intestine isoform,Peptide transporter 1,SLC15A1/PEPT1

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)