

— ABLYBIO, Help Your Research



# S100A5 Antibody

货号: **AYP6187**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	10kDa/12kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG
细胞定位	nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human S100A5.
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein has a Ca <sup>2+</sup> affinity 20- to 100-fold higher than the other S100 proteins studied under identical conditions. This protein also binds Zn <sup>2+</sup> and Cu <sup>2+</sup> , and Cu <sup>2+</sup> strongly which impairs the binding of Ca <sup>2+</sup> . This protein is expressed in very restricted regions of the adult brain.
基因ID	6276
基因名	S100A5
Swiss	P33763 ( <a href="https://www.uniprot.org/uniprotkb/P33763/entry">https://www.uniprot.org/uniprotkb/P33763/entry</a> )
别名	S100A5,S100D,S100A5 Antibody,Protein S-100D,S100 calcium-binding protein A5

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))