

— ABLYBIO, Help Your Research



RHOB Antibody

货号: **AYP6203**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	22kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HCT116,SK-BR-3,Daudi,Mouse brain,Rat lung
细胞定位	Cell membrane,Cleavage furrow,Late endosome membrane,Lipid-anchor,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human RHOB.
------	--

靶点信息

研究背景	RHOB (ras homolog family member B) is a protein-coding gene. GO annotations related to this gene include GDP binding and GTP binding. An important paralog of this gene is RND1.
------	--

基因ID	388
基因名	RHOB
Swiss	P62745
别名	RHOB,ARH6,ARHB,MST081,MSTP081,RHOH6,RHOB Antibody,Rho cDNA clone 6

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)