

— ABLYBIO, Help Your Research



# MRRF Antibody

货号: **AYP6206**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	24kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse testis,Mouse heart,Rat testis,Rat heart
细胞定位	mitochondrial matrix,mitochondrion
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human MRRF.
------	--

## 靶点信息

研究背景	This gene encodes a ribosome recycling factor, which is a component of the mitochondrial translational machinery. The encoded protein, along with mitochondrial elongation factor 2, functions in ribosomal recycling at the termination of mitochondrial translation by mediating the disassembly of ribosomes from messenger RNA. A pseudogene of this gene has been identified on chromosome X. [provided by RefSeq, Oct 2016]
基因ID	92399
基因名	MRRF
Swiss	Q96E11 ( <a href="https://www.uniprot.org/uniprotkb/Q96E11/entry">https://www.uniprot.org/uniprotkb/Q96E11/entry</a> )
别名	MRRF Antibody,MRRF,Ribosome-releasing factor,mitochondrial

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)