

— ABLYBIO, Help Your Research



SEPT6 Antibody

货号: **AYP6213**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF: 1:50 - 1:200
理论分子量	29kDa/48kDa/49kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Chromosome,Cleavage furrow,Cytoplasm,Midbody,centromere,cytoskeleton,kinetochore,spindle
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human SEPT6.
------	---

靶点信息

研究背景	This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis. One version of pediatric acute myeloid leukemia is the result of a reciprocal translocation between chromosomes 11 and X, with the breakpoint associated with the genes encoding the mixed-lineage leukemia and septin 2 proteins. This gene encodes four transcript variants encoding three distinct isoforms. An additional transcript variant has been identified, but its biological validity has not been determined.
基因ID	23157
基因名	SEPT6
Swiss	Q14141 (https://www.uniprot.org/uniprotkb/Q14141/entry)
别名	SEPT6,SEP2,SEPT2,septin-6,SEPT6 Antibody,KIAA0128

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)