

— ABLYBIO, Help Your Research



OXSR1 Antibody

货号: **AYP6218**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	58kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,293T,Jurkat
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human OXSR1.
------	---

靶点信息

研究背景	The product of this gene belongs to the Ser/Thr protein kinase family of proteins. It regulates downstream kinases in response to environmental stress, and may play a role in regulating the actin cytoskeleton.
------	---

基因ID	9943
基因名	OXSR1
Swiss	O95747 (https://www.uniprot.org/uniprotkb/O95747/entry)
别名	OXSR1,OSR1,OXSR1 Antibody,Oxidative stress-responsive 1 protein,KIAA1101

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)