

— ABLYBIO, Help Your Research



# RIOK3 Antibody

货号: **AYP6219**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	58kDa/59kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Mouse pancreas
细胞定位	cytosol
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human RIOK3.
------	---

## 靶点信息

研究背景	This gene was identified by the similarity of its product to the Aspergillus nidulans SUDD protein, an extra genic suppressor of the heat-sensitive bimD6 mutation that fails to attach properly to the spindle microtubules at a restrictive temperature. The specific function of this gene has not yet been determined.
------	--

基因ID	8780
基因名	RIOK3
Swiss	O14730 ( <a href="https://www.uniprot.org/uniprotkb/O14730/entry">https://www.uniprot.org/uniprotkb/O14730/entry</a> )
别名	RIOK3,SUDD,RIOK3 Antibody,RIO kinase 3,sudD homolog

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)