

— ABLYBIO, Help Your Research



SNX3 Antibody

货号: **AYP6235**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	13kDa/14kDa/16kDa/18kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,THP-1,SKOV3,Mouse kidney,Mouse spleen,Mouse brain,Rat kidney,Rat liver
细胞定位	Cytoplasmic vesicle,Early endosome,phagosome
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human SNX3.
------	--

靶点信息

研究背景	This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like most family members. This protein interacts with phosphatidylinositol-3-phosphate, and is involved in protein trafficking. A pseudogene of this gene is present on the sex chromosomes. Alternative splicing results in multiple transcript variants encoding distinct isoforms.
基因ID	8724
基因名	SNX3
Swiss	O60493 (https://www.uniprot.org/uniprotkb/O60493/entry)
别名	SNX3,Grd19,MCOPS8,SDP3,SNX3 Antibody,Protein SDP3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)