

— ABLYBIO, Help Your Research



CCT6A Antibody

货号: **AYP6249**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	53kDa/58kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human CCT6A.
------	---

靶点信息

研究背景	The protein encoded by this gene is a molecular chaperone that is a member of the chaperonin containing TCP1 complex (CCT), also known as the TCP1 ring complex (TRiC). This complex consists of two identical stacked rings, each containing eight different proteins. Unfolded polypeptides enter the central cavity of the complex and are folded in an ATP-dependent manner. The complex folds various proteins, including actin and tubulin. Alternate transcriptional splice variants of this gene, encoding different isoforms, have been characterized. In addition, several pseudogenes of this gene have been located.
基因ID	908
基因名	CCT6A
Swiss	P40227 (https://www.uniprot.org/uniprotkb/P40227/entry)
别名	CCT6A,CCT-zeta,CCT-zeta-1,CCT6,Cctz,HTR3,MoDP-2,TCP-1-zeta,TCP20,TCPZ,TTCP20,CCT6A Antibody,Acute morphine dependence-related protein 2,Chaperonin containing T-complex polypeptide 1 subunit 6A

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)