

— ABLYBIO, Help Your Research



UBA5 Antibody

货号: **AYP6269**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	38kDa/44kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,A-549,Jurkat,HeLa,Mouse brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human UBA5.
------	--

靶点信息

研究背景	This gene encodes a member of the E1-like ubiquitin-activating enzyme family. This protein activates ubiquitin-fold modifier 1, a ubiquitin-like post-translational modifier protein, via the formation of a high-energy thioester bond. Alternative splicing results in multiple transcript variants. A pseudogene of this gene has been identified on chromosome 1.
基因ID	79876
基因名	UBA5
Swiss	Q9GZZ9 (https://www.uniprot.org/uniprotkb/Q9GZZ9/entry)
别名	UBA5,EIEE44,SCAR24,THIFP1,UBE1DC1,UBA5 Antibody,UFM1-activating enzyme,Ubiquitin-activating enzyme E1 domain-containing protein 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)