

— ABLYBIO, Help Your Research



UBE2T Antibody

货号: **AYP6271**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	22kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Mouse thymus,Rat spleen,Rat heart
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human UBE2T.
------	---

靶点信息

研究背景	The protein encoded by this gene catalyzes the covalent attachment of ubiquitin to protein substrates. D effects in this gene have been associated with Fanconi anemia of complementation group T. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	29089
基因名	UBE2T
Swiss	Q9NPD8 (https://www.uniprot.org/uniprotkb/Q9NPD8/entry)
别名	UBE2T,FANCT,HSPC150,PIG50,UBE2T Antibody,Cell proliferation-inducing gene 50 protein,E2 ubiquitin-co njugating enzyme T,Ubiquitin carrier protein T,Ubiquitin-protein ligase T

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)