

— ABLYBIO, Help Your Research



WIPF1 Antibody

货号: **AYP6280**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	RAW264.7,Mouse spleen
细胞定位	actin cytoskeleton,cytoplasmic vesicle,cytosol
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human WIPF1.
------	---

靶点信息

研究背景	This gene encodes a protein that plays an important role in the organization of the actin cytoskeleton. The encoded protein binds to a region of Wiskott-Aldrich syndrome protein that is frequently mutated in Wiskott-Aldrich syndrome, an X-linked recessive disorder. Impairment of the interaction between these two proteins may contribute to the disease. Two transcript variants encoding the same protein have been identified for this gene.
基因ID	7456
基因名	WIPF1
Swiss	O43516 (https://www.uniprot.org/uniprotkb/O43516/entry)
别名	WIPF1,PRPL-2,WAS2,WASPIP,WIP,WIPF1 Antibody,Protein PRPL-2,Wiskott-Aldrich syndrome protein-interacting protein

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)