

— ABLYBIO, Help Your Research



SPON2 Antibody

货号: **AYP6317**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	35kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HCT116,HeLa,Rat brain
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human SPON2.
------	---

靶点信息

研究背景	Cell adhesion protein that promotes adhesion and outgrowth of hippocampal embryonic neurons. Binds directly to bacteria and their components and functions as an opsonin for macrophage phagocytosis of bacteria. Essential in the initiation of the innate immune response and represents a unique pattern-recognition molecule in the ECM for microbial pathogens (By similarity. Binds bacterial lipopolysaccharide (LPS.
基因ID	10417
基因名	SPON2
Swiss	Q9BUD6 (https://www.uniprot.org/uniprotkb/Q9BUD6/entry)
别名	SPON2,DIL-1,DIL1,M-SPONDIN,MINDIN,spondin-2,SPON2 Antibody,Differentially expressed in cancerous and non-cancerous lung cells 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)