

— ABLYBIO, Help Your Research



IGF2BP2 Antibody

货号: **AYP6338**

产品信息

| | |
|-------|---|
| 反应 | Human |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IF ELISA |
| 推荐浓度 | WB: 1:500 - 1:2000 IF: 1:50 - 1:200 |
| 理论分子量 | 54kDa/58kDa/59kDa/61kDa/66kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | HT-29,BxPC-3,HeLa,Jurkat,NIH/3T3 |
| 细胞定位 | Cytoplasm,Nucleus |
| 纯化 | Affinity purification |

抗原信息

| | |
|------|---|
| 抗原信息 | Synthesized peptide derived from Human IGF2BP2. |
|------|---|

靶点信息

| | |
|-------|--|
| 研究背景 | This gene encodes a member of the IGF-II mRNA-binding protein (IMP) family. The protein encoded by this gene contains four KH domains and two RRM domains. It functions by binding to the 5' UTR of the insulin-like growth factor 2 (IGF2) mRNA and regulating IGF2 translation. Alternative promoter usage and alternative splicing result in multiple variants encoding different isoforms. |
| 基因ID | 10644 |
| 基因名 | IGF2BP2 |
| Swiss | Q9Y6M1 (https://www.uniprot.org/uniprotkb/Q9Y6M1/entry) |
| 别名 | IGF2BP2,IMP-2,IMP2,VICKZ2,IGF2BP2 Antibody,Hepatocellular carcinoma autoantigen p62,IGF-II mRNA-binding protein 2,VICKZ family member 2 |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)