

— ABLYBIO, Help Your Research



IF16 Antibody

货号: **AYP6341**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	12kDa/13kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human IF16.
------	--

靶点信息

研究背景	This gene was first identified as one of the many genes induced by interferon. The encoded protein may play a critical role in the regulation of apoptosis. A minisatellite that consists of 26 repeats of a 12 nucleotide repeating element resembling the mammalian splice donor consensus sequence begins near the end of the second exon. Alternatively spliced transcript variants that encode different isoforms by using the two downstream repeat units as splice donor sites have been described.
基因ID	2537
基因名	IFI6
Swiss	P09912 (https://www.uniprot.org/uniprotkb/P09912/entry)
别名	IFI6,6-16,FAM14C,G1P3,IFI-6-16,IFI616,IFI6 Antibody,Interferon-induced protein 6-16

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)