

— ABLYBIO, Help Your Research



# CYP4X1 Antibody

货号: **AYP6375**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney
细胞定位	Endoplasmic reticulum membrane
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human CYP4X1.
------	--

## 靶点信息

研究背景	This gene encodes a member of the cytochrome P450 superfamily of enzymes and is located within a cluster of genes belonging to this superfamily on chromosome 1. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. The expression pattern of a similar rat protein suggests that this protein may be involved in neurovascular function in the brain. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]
基因ID	260293
基因名	CYP4X1
Swiss	Q8N118 ( <a href="https://www.uniprot.org/uniprotkb/Q8N118/entry">https://www.uniprot.org/uniprotkb/Q8N118/entry</a> )
别名	CYP4X1, CYP4X1 Antibody

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))