

— ABLYBIO, Help Your Research



# MERTK Antibody

货号: **AYP6400**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	110kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse spleen
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human MERTK.
------	---

## 靶点信息

研究背景	This gene is a member of the MER/AXL/TYRO3 receptor kinase family and encodes a transmembrane protein with two fibronectin type-III domains, two Ig-like C2-type (immunoglobulin-like) domains, and one tyrosine kinase domain. Mutations in this gene have been associated with disruption of the retinal pigment epithelium (RPE) phagocytosis pathway and onset of autosomal recessive retinitis pigmentosa (RP).
基因ID	10461
基因名	MERTK
Swiss	Q12866 ( <a href="https://www.uniprot.org/uniprotkb/Q12866/entry">https://www.uniprot.org/uniprotkb/Q12866/entry</a> )
别名	MERTK,MER,RP38,Tyro12,c-Eyk,c-mer,MERTK Antibody,Proto-oncogene c-Mer,Receptor tyrosine kinase MerTK

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))