

— ABLYBIO, Help Your Research



# APOE Antibody

货号: **AYP6408**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	36kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Mouse brain
细胞定位	Secreted
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human APOE.
------	--

## 靶点信息

研究背景	This gene encodes a member of the apolipoprotein A1/A4/E family of proteins. This protein is involved in the transport of lipoproteins in the blood. It binds to a specific liver and peripheral cell receptor, and is essential for the normal catabolism of triglyceride-rich lipoprotein constituents. Homozygous knockout mice for this gene accumulate high levels of cholesterol in the blood and develop atherosclerosis. Different alleles of this gene have been associated with either increased risk or a protective effect for Alzheimer's disease in human patients. This gene maps to chromosome 7 in a cluster with the related apolipoprotein C1, C2 and C4 genes. [provided by RefSeq, Apr 2015]
基因ID	11816
基因名	APOE
Swiss	P02649 ( <a href="https://www.uniprot.org/uniprotkb/P02649/entry">https://www.uniprot.org/uniprotkb/P02649/entry</a> )
别名	APOE,AD2,APO-E,ApoE4,LDLCQ5,LPG,AI255918,APOE Antibody

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)