

— ABLYBIO, Help Your Research



FCGR2A/CD32 Antibody

货号: **AYP6410**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	34kDa/35kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,Mouse spleen
细胞定位	Cell membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human FCGR2A/CD32.
------	---

靶点信息

研究背景	This gene encodes one member of a family of immunoglobulin Fc receptor genes found on the surface of many immune response cells. The protein encoded by this gene is a cell surface receptor found on phagocytic cells such as macrophages and neutrophils, and is involved in the process of phagocytosis and clearing of immune complexes. Alternative splicing results in multiple transcript variants.
基因ID	2212
基因名	FCGR2A
Swiss	P12318 (https://www.uniprot.org/uniprotkb/P12318/entry)
别名	FCGR2A,CD32,CD32A,CDw32,FCG2,FCGR2,FCGR2A1,FcGR,IGFR2,FCGR2A/CD32 Antibody,Fc-gamma RII-a,FCGR2A/CD32

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)