

— ABLYBIO, Help Your Research



# SLC44A1 Antibody

货号: **AYP6414**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	73kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,LO2,mouse liver,mouse brain,mouse kidney,rat liver
细胞定位	Cell membrane,Mitochondrion outer membrane,Multi-pass membrane protein
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human SLC44A1.
------	---

## 靶点信息

研究背景	Choline transporter. Involved in membrane synthesis and myelin production.
------	--

基因ID	23446
基因名	SLC44A1
Swiss	Q8WWI5 ( <a href="https://www.uniprot.org/uniprotkb/Q8WWI5/entry">https://www.uniprot.org/uniprotkb/Q8WWI5/entry</a> )
别名	SLC44A1,CD92,CDW92,CHTL1,CTL1,SLC44A1 Antibody,Solute carrier family 44 member 1

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)