

— ABLYBIO, Help Your Research



# CLM-1 Antibody

货号: **AYP6421**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	18kDa/21kDa/26kDa/32kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,U-937,Rat brain
细胞定位	Cell membrane,Single-pass type I membrane protein
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human CLM-1.
------	---

## 靶点信息

研究背景	This gene encodes a member of the CD300 protein family. Members of this family are cell surface glycoproteins with a single IgV-like extracellular domain, and are involved in the regulation of immune response. The encoded protein is an inhibitory receptor. Alternative splicing results in multiple transcript variants.
------	--

基因ID	146722
基因名	CD300LF
Swiss	Q8TDQ1
别名	CD300LF,CD300f,CLM-1,CLM1,IREM-1,IREM1,IgSF13,LMIR3,NKIR,CLM-1 Antibody,CD300 antigen-like family member F,Immune receptor expressed on myeloid cells 1,Immunoglobulin superfamily member 13,NK inhibitory receptor

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))