

# DSG1 Antibody

货号: **AYP6428**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	42kDa/113kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,HeLa,293T,BT-474
细胞定位	Cell junction,Cell membrane,Single-pass type I membrane protein,desmosome
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human DSG1.
------	--

## 靶点信息

研究背景	This gene encodes a member of the desmoglein protein subfamily. Desmogleins, along with desmocollins, are cadherin-like transmembrane glycoproteins that are major components of the desmosome. Desmosomes are cell-cell junctions that help resist shearing forces and are found in high concentrations in cells subject to mechanical stress. This gene is found in a cluster with other desmoglein family members on chromosome 18. The encoded protein has been identified as a target of auto-antibodies in the autoimmune skin blistering disease pemphigus foliaceus. Disruption of this gene has also been associated with the skin diseases palmoplantar keratoderma and erythroderma.
基因ID	1828
基因名	DSG1

<b>Swiss</b>	Q02413
别名	DSG1;CDHF4;DG1;DSG;EPKHE;EPKHIA;PPKS1;SPPK1

产品验证

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)