

— ABLYBIO, Help Your Research



# FTH1 Antibody

货号: **AYP6432**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	21kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	cytoplasm,cytosol,extracellular exosome,extracellular region,nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human FTH1.
------	--

## 靶点信息

研究背景	This gene encodes the heavy subunit of ferritin, the major intracellular iron storage protein in prokaryotes and eukaryotes. It is composed of 24 subunits of the heavy and light ferritin chains. Variation in ferritin subunit composition may affect the rates of iron uptake and release in different tissues. A major function of ferritin is the storage of iron in a soluble and nontoxic state. Defects in ferritin proteins are associated with several neurodegenerative diseases. This gene has multiple pseudogenes. Several alternatively spliced transcript variants have been observed, but their biological validity has not been determined.
基因ID	2495
基因名	FTH1
Swiss	P02794
别名	FHC,FTH,FTHL6,HFE5,PIG15,PLIF,Ferritin,FTH1,FTH1 Antibody,Cell proliferation-inducing gene 15 protein

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)