

— ABLYBIO, Help Your Research



MAD2L1 Antibody

货号: **AYP6441**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	10kDa/23kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Daudi cells, HeLa cells, HEK-293 cells, rat liver tissue
细胞定位	Chromosome,Cytoplasm,Nucleus,centromere,cytoskeleton,kinetochore,spindle pole
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human MAD2L1.
------	--

靶点信息

研究背景	MAD2L1 is a component of the mitotic spindle assembly checkpoint that prevents the onset of anaphase until all chromosomes are properly aligned at the metaphase plate. MAD2L1 is related to the MAD2L2 gene located on chromosome 1. A MAD2 pseudogene has been mapped to chromosome 14.
------	---

基因ID	4085
基因名	MAD2L1
Swiss	Q13257 (https://www.uniprot.org/uniprotkb/Q13257/entry)
别名	MAD2L1,HSMAD2,MAD2,MAD2L1 Antibody,Mitotic arrest deficient 2-like protein 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)