

— ABLYBIO, Help Your Research



EDG3 Antibody

货号: **AYP6462**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	42kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Mouse kidney
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human EDG3.
------	--

靶点信息

研究背景	This gene encodes a member of the EDG family of receptors, which are G protein-coupled receptors. This protein has been identified as a functional receptor for sphingosine 1-phosphate and likely contributes to the regulation of angiogenesis and vascular endothelial cell function.
基因ID	1903
基因名	S1PR3
Swiss	Q99500 (https://www.uniprot.org/uniprotkb/Q99500/entry)
别名	S1PR3,EDG-3,EDG3,LPB3,S1P3,EDG3 Antibody,Endothelial differentiation G-protein coupled receptor 3,Sphingosine 1-phosphate receptor Edg-3,C9orf108

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)