

— ABLYBIO, Help Your Research



Cleaved Cathepsin L2 (Leu114) Rabbit Antibody

货号: **AYP6535**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	37kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,Mouse liver,Rat liver
细胞定位	Lysosome
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Cleaved Cathepsin L2 (Leu114).
------	---

靶点信息

研究背景	The protein encoded by this gene, a member of the peptidase C1 family, is a lysosomal cysteine proteinase that may play an important role in corneal physiology. This gene is expressed in colorectal and breast carcinomas but not in normal colon, mammary gland, or peritumoral tissues, suggesting a possible role for this gene in tumor processes. Alternatively spliced variants, encoding the same protein, have been identified.
基因ID	1515
基因名	CTSV
Swiss	O60911
别名	CTSV,CATL2,CTSL2,CTSU,Cleaved Cathepsin L2 (Leu114) Rabbit Antibody,Cathepsin U,Cathepsin V,Cleaved Cathepsin L2 (Leu114)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)