

— ABLYBIO, Help Your Research



# Cleaved C1s HC (Arg437) Antibody

货号: **AYP6547**

## 产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	76kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	blood microparticle,complement component C1 complex,extracellular region,extracellular space
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human Cleaved C1s HC (Arg437).
------	---

## 靶点信息

研究背景	This gene encodes a serine protease, which is a major constituent of the human complement subcomponent C1. C1s associates with two other complement components C1r and C1q in order to yield the first component of the serum complement system. Defects in this gene are the cause of selective C1s deficiency.
基因ID	716
基因名	C1S
Swiss	P09871
别名	C1S,EDSPD2,Cleaved C1s HC (Arg437) Antibody,C1 esterase,Complement component 1 subcomponent s, Complement C1s subcomponent chain A,Complement C1s subcomponent chain B,Cleaved C1s HC (Arg437)

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)