

— ABLYBIO, Help Your Research



Iba 1 Mouse mAb

货号: **AYS1001**

产品信息

反应	Human,Mouse,Rat
宿主	Mouse
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB 1:500-2000;IHC 1:100-700 IF 1:100-500;:
理论分子量	17kD(Calculated).
实测分子量	17kDa;
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated.
阳性对照	
细胞定位	Cytoplasm>Cytoskeleton. Cell projection>Ruffle membrane>Peripheral membrane protein>Cytoplasmic side. Cell projection>Phagocytic cup.
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from part region of human protein AA range: 70-147
------	--

靶点信息

研究背景	Actin-binding protein that enhances membrane ruffling and RAC activation. Enhances the actin-bundling activity of LCP1. Binds calcium. Plays a role in RAC signaling and in phagocytosis. May play a role in macrophage activation and function. Promotes the proliferation of vascular smooth muscle cells and of T-lymphocytes. Enhances lymphocyte migration. Plays a role in vascular inflammation
基因ID	199;
基因名	AIF1 G1 IBA1
Swiss	P55008
别名	AIF 1; AIF-1; Aif1; AIF1 protein; AIF1_HUMAN; Allograft inflammatory factor 1; Allograft inflammatory factor 1 splice variant G; allograft inflammatory factor-1 splice variant Hara-1; balloon angioplasty responsive transcription; BART 1; G1; G1 putative splice variant of allograft inflammatory factor 1; IBA 1; IBA1; interferon gamma responsive transcript; Interferon responsive transcript 1; interferon responsive transcript factor 1; Ionized calcium binding adapter molecule 1; Ionized calcium-binding adapter molecule 1; ionized calcium-binding adapter molecule; IRT 1; IRT1; Microglia response factor; MRF1; Protein g1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)