

— ABLYBIO, Help Your Research



# NPNT Rabbit pAb

货号: **AYT1005**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB 1:500-2000:</b>
理论分子量	62 kDa
实测分子量	62-70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted>Extracellular space>Extracellular matrix.
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from human NPNT AA range: 470-538
------	---

## 靶点信息

研究背景	Nephronectin also known as NPNT is an extracellular matrix protein that functions mechanically as a ligand for integrin $\alpha 8\beta 1$ and integrin $\alpha 4\beta 1$ . Nephronectin has a mass of approximately 70 kDa. It is mainly expressed in the kidneys specifically within the glomerular basement membrane and mesangium. Additionally it can be found in tissues such as the lungs heart and bone. Nephronectin plays a role in promoting cell adhesion and migration particularly in epithelial cells and osteoblasts.
基因ID	255743
基因名	NPNT EGFL6L
Swiss	Q6UXI9 ( <a href="https://www.uniprot.org/uniprotkb/Q6UXI9/entry">https://www.uniprot.org/uniprotkb/Q6UXI9/entry</a> )
别名	展开/折叠 EGFL6L; Nephronectin; NPNT; NPNT_HUMAN; POEM; Preosteoblast EGF like repeat protein with MAM domain; Preosteoblast EGF-like repeat protein with MAM domain; Protein EGFL6 like; Protein EGFL6-like;

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)