



RPS16 Antibody

货号: **AYT1034**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ICC/IF
推荐浓度	WB 1:500-2000: IF, 1: 50— 1: 300:
理论分子量	16KD
实测分子量	16KD
形式	Liquid
保存条件	Store at -20℃. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.35.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	RPS16 also known as Ribosomal Protein S16 is a component of the 40S ribosomal subunit with a molecular mass of approximately 16.5 kDa. It plays a role in protein synthesis by joining together with other ribosomal proteins and rRNA to form the ribosome's small subunit. RPS16 is expressed in various tissues throughout the human body reflecting its essential role in maintaining protein synthesis machinery.
------	--

基因ID	6217
基因名	RPS16
Swiss	P62249 (https://www.uniprot.org/uniprotkb/P62249/entry)
别名	40S ribosomal protein S16; Ribosomal protein S16; RPS 16; S16;

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)