



IL-1β Mouse mAb

货号: **AYT1036**

产品信息

反应	Human,Rat,Mouse
宿主	Mouse
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB 1:1000-1:5000: IHC 1:200-1000: IF 1:100-500:
理论分子量	31kD
实测分子量	35kD,31kD,19kD
形式	Liquid
保存条件	Store at -20℃. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.135.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	亲和纯化


抗原信息

抗原信息	
------	--


靶点信息

研究背景	interleukin 1 beta (IL1B) Homo sapiens The protein encoded by this gene is a member of the interleukin 1 cytokine family. This cytokine is produced by activated macrophages as a proprotein , which is proteolytic ally processed to its active form by caspase 1 (CASP1/ICE) . This cytokine is an important mediator of the inflammatory response , and is involved in a variety of cellular activities , including cell proliferation , differentiation , and apoptosis. The induction of cyclooxygenase-2 (PTGS2/COX2) by this cytokine in the central nervous system (CNS) is found to contribute to inflammatory pain hypersensitivity. This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2
基因ID	3553;
基因名	IL1B IL1F2
Swiss	P01584 (https://www.uniprot.org/uniprotkb/P01584/entry)
别名	IL 1; IL 1 beta; IL-1 beta; IL1 BETA; IL1B; IL1B_HUMAN; IL1F2; Interleukin 1 beta; Interleukin-1 beta; OAF; OTTHUMP00000162031; Preinterleukin 1 beta; Pro interleukin 1 beta;

产品验证



Western blot analysis of extracts from mouse lung(LPS),RAW 264.7(LPS) cells using IL-1β at 1:1000.



Immunohistochemical analysis of paraffin-embedded Mouse-heart tissue., IL-1β Mouse mAb

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)