



KRT35 Rabbit pAb

货号: **AYT1039**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	50 kDa
实测分子量	50 kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated.
阳性对照	
细胞定位	extracellular space ,intermediate filament ,extracellular exosome ,
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from human protein . at AA range: 305-411
------	---

靶点信息

研究背景	The protein encoded by this gene is a member of the keratin gene family. This type I hair keratin is an aci dic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. [provided b y RefSeq
------	---

基因ID	3886;
基因名	KRT35 HHA5 HKA5 KRTHA5
Swiss	Q92764 (https://www.uniprot.org/uniprotkb/Q92764/entry)
别名	Hair keratin; Hair keratin, type I Ha5; HHA5; HKA5; K35; Keratin; Keratin, type I cuticular Ha5; Keratin-35; KRT35; KRT35_HUMAN; KRTHA5; type I cuticular Ha5; type I Ha5;

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)