



Wnt5b Rabbit mAb

货号: **AYT1043**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:100-1:700 IF 1:100-1:700: ELISA 1:5000-1:20000:
理论分子量	40kDa
实测分子量	45KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse ovary
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 120-359 of human WNT 5B (NP_116031.1).
------	--


靶点信息

研究背景	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. The se proteins have been implicated in oncogenesis and in several developmental processes, including regul ation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 94% and 80% amino acid identity to the mouse Wnt5b protein and the hu man WNT5A protein, respectively. Alternative splicing of this gene generates 2 transcript variants.
基因ID	81029
基因名	WNT5B
Swiss	Q9H1J7 (https://www.uniprot.org/uniprotkb/Q9H1J7/entry)
别名	WNT5B,WNT5B Rabbit pAb

产品验证



Western blot analysis of various lysates using Wnt5b at 1:1000 dilution.



Human tonsil was stained with anti-Wnt5b Rabbit antibody, 1: 300

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)