



ATG101 Rabbit pAb

货号: **AYT1047**

产品信息

反应	Human,Mouse,Rat
宿主	Mouse
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF
推荐浓度	WB 1:500-2000: IF/ICC 1:100-1:500:
理论分子量	24kD
实测分子量	25kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.831.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	ATG101 also known as Autophagy Related Protein 101 is a critical component of the cellular process known as autophagy. This protein has an approximate mass of 289 amino acids making it a modest-sized protein. ATG101 is expressed in various cell types especially in those with high autophagic activity. It mainly resides in the cytoplasm where it participates in the initiation and regulation of autophagosome formation.
基因ID	60673
基因名	ATG101 C12orf44 PP894
Swiss	Q9BSB4 (https://www.uniprot.org/uniprotkb/Q9BSB4/entry)
别名	ATG101; Atg13-interacting protein; ATGA1_HUMAN; Autophagy-related protein 101; C12orf44; Chromosome 12 open reading frame 44; FLJ11773; OTTHUMP00000241687; OTTHUMP00000241688; OTTHUMP00000241689;

产品验证

Western blot analysis of extracts from HUVEC tissue using ATG101 at 1:500.



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)