



ZEB1 Rabbit mAb

货号: **AYT1049**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB: WB 1:1000-1:5000 IHC 1:100-1:500: IF 1:200-1:800: ELISA 1:5000-1:20000:
理论分子量	116kDa/122kDa/124kDa
实测分子量	200KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,Rat heart
细胞定位	cytosol,nucleoplasm,nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant protein of human ZEB1
------	-----------------------------------


靶点信息

研究背景	This gene encodes a zinc finger transcription factor. The encoded protein likely plays a role in transcriptional repression of interleukin 2. Mutations in this gene have been associated with posterior polymorphous corneal dystrophy-3 and late-onset Fuchs endothelial corneal dystrophy. Alternatively spliced transcript variants encoding different isoforms have been described.
基因ID	6935
基因名	ZEB1
Swiss	P37275 (https://www.uniprot.org/uniprotkb/P37275/entry)
别名	ZEB1,AREB6,BZP,DELTAEF1,FECD6,NIL2A,PPCD3,TCF8,ZFHEP,ZFHX1A,ZEB1 Rabbit pAb,NIL-2-A zinc finger protein,Negative regulator of IL2,Transcription factor 8

产品验证



Western blot analysis of ZEB1 expressed in HeLa,Jurkat,Rat heart using ZEB1 Rabbit at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.



Human brain was stained with anti-ZEB1 Rabbit antibody, 1: 100

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)