



DHRS1 Rabbit pAb

货号: **AYT1061**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	
理论分子量	34 kDa
实测分子量	34 kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm, Nucleus, Cell junction, tight junction, Cell membrane
纯化	亲和纯化

抗原信息


抗原信息	Synthesized peptide derived from the Internal region of human DHRS1
------	---

靶点信息

研究背景	DHRS1(Dehydrogenase/reductase SDR family member 1) is also named as LJ25430, MGC20204, DKFZp586I0523, SDR19C1, FLJ14250 and belongs to the short-chain dehydrogenases/reductases (SDR) family. It is involved in the metabolism of a large variety of compounds, including steroid hormones, prostaglandins, retinoids, lipids and xenobiotics.Northern blot analysis detected DHRS1express in heart highest and in liv er lowes
------	--

基因ID	115817;
基因名	DHRS1
Swiss	Q96LJ7 (https://www.uniprot.org/uniprotkb/Q96LJ7/entry)
别名	展开/折叠 1110029G07Rik; AW112170; D14Ertd484e; dehydrogenase/reductase (SDR family) member 1; Dehydrogenase/reductase SDR family member 1; DHRS; Dhrrs1; DHRS1_HUMAN; FLJ25430; MGC20204; SDR19C1;

产品验证



Western blot detection of DHRS1 in C6 cell lysates using DHRS1 antibody(1:500 diluted).

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)