



GFAP Mouse mAb

货号: **AYT1063**

产品信息

反应	Human,Mouse,Rat
宿主	Mouse
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB 1:1000-5000; : IHC 1:50-300: IF 1:100-1:300:
理论分子量	50kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20℃. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.302.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm . Associated with intermediate filaments.
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	GFAP (Glial fibrillary acidic protein) is a type III intermediate filament (IF) protein specific to the central nervous system (CNS). GFAP is one of the main components of the intermediate filament network in astrocytes and has been proposed as playing a role in cell migration, cell motility, maintaining mechanical strength, and in mitosis.
基因ID	2670;
基因名	GFAP
Swiss	P14136 (https://www.uniprot.org/uniprotkb/P14136/entry)
别名	wu:fb34h11; ALXDRD; cb345; etID36982.3; FLJ42474; FLJ45472; GFAP; GFAP_HUMAN; gfapl; Glial fibrillary acidic protein; Intermediate filament protein; wu:fk42c12; xx:af506734; zgc:110485;

产品验证



Human cerebrum tissue was stained with GFAP Antibody,1:100

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)