



NeuN Mouse mAb

货号: **AYT1064**

产品信息

反应	Human,Mouse
宿主	Mouse
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB 1:500-2000: IHC 1:100-1000,; .IF 1:100-500:
理论分子量	34kD
实测分子量	45-55kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.1029.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasmic , Nuclear
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	NeuN also known as Fox-3 is a protein often utilized as a neuronal marker due to its expression predominantly in neurons. The molecular weight of NeuN is approximately 46-48 kDa. It is located in the nucleus and in the post-mitotic neurons of the central nervous system. Research frequently uses anti-NeuN antibodies for identifying neuronal cells in various samples capitalizing on its reliable expression profile across diverse neuron types. Common techniques for detecting NeuN include Western blot and immunohistochemistry (IHC) with NeuN staining being a definitive method for confirming neuronal identity. Biological function summary
基因ID	146713;
基因名	RBFOX3
Swiss	A6NFN3 (https://www.uniprot.org/uniprotkb/A6NFN3/entry)
别名	FLJ56884; FLJ58356; Fox-1 homolog C; fox1 homolog C; Fox3; FOX3NeuN; hexaribonucleotide binding protein 3; HRNBP3; NEUN; neuronal nuclei; Rbfox3; RFOX3_HUMAN; RNA binding protein fox-1 homolog 3; RNA binding protein, fox 1 homolog (C. elegans) 3;

产品验证



Mouse brain was stained with anti-nuclear NeuN Mouse antibody,1:200

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)