



TMED1 Rabbit pAb

货号: **AYT1067**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ICC/IF
推荐浓度	WB 1:500-2000;; IF: 1: 50-1: 200:
理论分子量	25KD
实测分子量	30 kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell membrane ; Single-pass type I membrane protein
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Transmembrane p24 trafficking protein 1 (TMED1) belongs to the TMED protein family, also known as the p24 family, and plays a part in protein vesicular trafficking. As the critical regulator of protein trafficking along the secretory pathway, the TMED proteins influence the composition, structure, and function of the endoplasmic reticulum and the Golgi apparatus. Disruption of the secretory pathway could affect cell differentiation and survival, ultimately resulting in disease
基因ID	11018;
基因名	TMED1 IL1RL1L IL1RL1LG
Swiss	Q13445 (https://www.uniprot.org/uniprotkb/Q13445/entry)
别名	IL1RL1L, IL1RL1LG, Interleukin-1 receptor-like 1 ligand, p24 family protein gamma-1, p24gamma1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)