

E-cadherin Mouse Monoclonal Antibody

货号: **AYT1069**

产品信息

反应	Human,Mouse,Rat
宿主	Mouse
克隆性	Monoclonal
预测反应	
应用	ICC IHC ICC/IF FC
推荐浓度	WB : 1/500 - 1/2000 IHC : 1/100 - 1/1000 IF/ICC : 1/100 - 1/1000 FC : 1/200 - 1/400
理论分子量	97kD(Calculated).
实测分子量	80-120kDa(cleavages),135 kDa(Observed)
形式	
保存条件	
偶联物	
阳性对照	
细胞定位	
纯化	

抗原信息


抗原信息	
------	--

靶点信息


研究背景	
------	--

基因ID	CDH1(999)
基因名	CDH1
Swiss	P12830 (https://www.uniprot.org/uniprotkb/P12830/entry)
别名	Arc 1; CADH1_HUMAN; Cadherin 1; cadherin 1 type 1 E-cadherin; Cadherin1; CAM 120/80; CD 324; CD324; CD324 antigen; cdh1; CDHE; E-Cad/CTF3; E-cadherin; ECAD; Epithelial cadherin; epithelial calcium dependent adhesion protein; LCAM; Liver cell adhesion molecule; UVO; Uvomorulin;


产品验证



Western blot analysis using CDH1 mouse mAb against PC-3 (1), A431 (2), MCF-7 (3), HT-29 (4), HePG2 (5), C6 (6), A549 (7), NIH/3T3 (8), and C2C12 (9) cell lysate.



Immunohistochemical analysis of paraffin embedded human stomach cancer tissues using CDH1 mouse mAb with DAB staining.1:200



Immunofluorescence analysis of Hela cells using CDH1 mouse mAb ,1:200

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)