



SUN5 Rabbit pAb

货号: **AYT1076**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ICC/IF
推荐浓度	WB 1:500-1:2000: IF/ICC 1:100-1:500:
理论分子量	41kD
实测分子量	43kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Mitochondrion inner membrane>Multi-pass membrane protein.
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from part region of human protein
------	---

靶点信息

研究背景	SPAG4L,also named as SUN domain-containing protein 5 or Testis and spermatogenesis-related gene 4 protein, is a 379 amino acid protein, which contains 1 SUN domain. SPAG4L is widely expressed and localizes in the nucleus inner membrane. SPAG4L is restricted to the apical nuclear region of round spermatids that face the acrosomic vesicle. SPAG4L is involved in linkage of the acrosomic vesicle to the spermatid nucleus, and in acrosome biogenesis.
基因ID	140732;
基因名	SUN5 SPAG4L TSARG4
Swiss	Q8TC36 (https://www.uniprot.org/uniprotkb/Q8TC36/entry)
别名	dj726C3.1; MGC33594; SPAG4-like protein; Sperm associated antigen 4 like protein; SUN domain-containing protein 5; SUN5; SUN5 Sad1 and UNC84 domain containing 5; testis and spermatogenesis related gene 4; Testis and spermatogenesis related gene 4 protein; TSARG4;

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)