



IFT20 Rabbit pAb

货号: **AYT1101**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB 1:500-1:2000, : IF: 1:50-1:200: IHC 1:50-1:200:
理论分子量	15 kDa
实测分子量	15-18 kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	
阳性对照	
细胞定位	Golgi apparatus>cis-Golgi network. Cytoplasm>Cytoskeleton>Microtubule organizing center>Centrosome>Centriole. Cytoplasm>Cytoskeleton>Cilium basal body. Cell projection>Cilium. Cytoplasm>Cytoskeleton. Golgi apparatus.
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from human IFT20 AA
------	---

靶点信息

研究背景	Ifi204 also known as Interferon-activated gene 204 is a protein of the mouse involved in the immune response. It has a molecular mass of approximately 74 kilodaltons. This protein is expressed mainly in activated macrophages and neutrophils and contributes to the innate immune mechanism. It detects cytosolic DNA and activates the STING pathway to induce type I interferon production. It contains a HIN-200 domain that helps in recognizing specific DNA sequences.
基因ID	
基因名	IFT20
Swiss	Q8IY31 (https://www.uniprot.org/uniprotkb/Q8IY31/entry)
别名	hIFT20; Intraflagellar transport 20 homolog (Chlamydomonas); Intraflagellar transport protein 20 like; Intraflagellar transport protein IFT20

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)