



SAA3 Rabbit pAb

货号: **AYT1104**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	13KD
实测分子量	13KD-15KD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	无明确组织表达
细胞定位	cytoplasmic
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Mouse SAA3
------	---

靶点信息

研究背景	Serum Amyloid A3 (SAA3) is a member of the serum amyloid A family of proteins and is known for its acute-phase response activity. This protein has an approximate mass of 14 kDa. SAA3 expression is mainly found in extrahepatic tissues in response to pro-inflammatory stimuli. It increases quickly during inflammation and tissue injury serving as a sensitive marker for inflammatory diseases. SAA3's function in the immune response makes it an important protein for researchers studying inflammation and related conditions.
基因ID	20210;
基因名	Saa3
Swiss	P04918 (https://www.uniprot.org/uniprotkb/P04918/entry)
别名	Serum amyloid A-3 protein, SAA3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)