



CD80 (YD16115) Rabbit mAb

货号: **AYT1105**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	
理论分子量	65kD
实测分子量	60-70kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Cell membrane
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息


研究背景	The protein encoded by this gene is a membrane receptor that is activated by the binding of CD28 or CTL A-4. The activated protein induces T-cell proliferation and cytokine production. This protein can act as a receptor for adenovirus subgroup B and may play a role in lupus neuropathy.
------	--

基因ID	941
基因名	CD80
Swiss	P33681 (https://www.uniprot.org/uniprotkb/P33681/entry)
别名	CD80 (YD16115),CD80 (YD16115) Rabbit mAb,CD80,Activation B7-1 antigen,BB1,CTLA-4 counter-receptor B7.1,CD28LG,CD28LG1,LAB7


产品验证



Western blot analysis of extracts from Daudi Cell using CD80 at 1:1000.



Human Lymph node was stained with anti-CD80 rabbit antibody, 1: 100



Immunofluorescence analysis of RAW264.7 (LPS stimulation), CD80 (YD16115) Rabbit mAb, 1: 100

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)