



RNAS4 Rabbit pAb

货号: **AYT1108**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB 1:500-2000: IHC 1:50-1:200:
理论分子量	16kD
实测分子量	17-20 kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted.
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from human protein
------	--

靶点信息

研究背景	Ribonuclease 4 also known as RNase 4 is a member of the ribonuclease A superfamily. This protein has a n approximate mass of 16 kDa. RNase 4 functions as an endonuclease catalyzing the cleavage of RNA on the 3' side of pyrimidines like uridine and cytidine. It shows high expression levels in the pancreas and liver and is also found in other tissues such as the spleen and heart
基因ID	6038;
基因名	RNASE4 RNS4
Swiss	P34096 (https://www.uniprot.org/uniprotkb/P34096/entry)
别名	Ribonuclease 4, RNase 4, RNASE4, RNS4, ribonuclease, RNase A family, 4

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)