



Gas1 Rabbit pAb

货号: **AYT1113**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB 1:500-1:1000, IF/ICC 1:100-1:500, IHC 1:50-1:200:
理论分子量	36 KD
实测分子量	36 KD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	无明确组织表达
细胞定位	Cell membrane
纯化	Affinity purification

抗原信息

抗原信息	The antiserum was produced against synthesized peptide derived from human GAS1
------	--

靶点信息

研究背景	The growth arrest-specific 1 protein or Gas1 plays an inhibitory role in cell cycle progression by interacting with various growth factors and receptors. Gas1 often termed as Gas-1 is a 38 kDa protein. It expresses prominently in quiescent cells particularly in the nervous and embryonic tissues. This expression helps control the proliferation of cells which is critical for maintaining normal cellular functions and structure.
------	--

基因ID	2619;
基因名	Growth arrest-specific protein 1
Swiss	P54826 (https://www.uniprot.org/uniprotkb/P54826/entry)
别名	GAS-1; Gas1; GAS1_HUMAN; Growth arrest-specific protein 1; Growth-arrest-specific protein 1 precursor;

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)