



MBL2 Antibody

货号: **AYT1114**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB: 1:1000-1:5000 IHC 1:100-1:600: IF 1:100-1:700:
理论分子量	26kDa
实测分子量	26kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse serum
细胞定位	Secreted
纯化	Affinity purification


抗原信息

抗原信息	Synthesized peptide derived from Human MBL2.
------	--


靶点信息

研究背景	This gene encodes the soluble mannose-binding lectin or mannose-binding protein found in serum. The p rotein encoded belongs to the collectin family and is an important element in the innate immune system. The protein recognizes mannose and N-acetylglucosamine on many microorganisms, and is capable of ac tivating the classical complement pathway. Deficiencies of this gene have been associated with susceptib ility to autoimmune and infectious diseases.
基因ID	4153
基因名	MBL2 COLEC1 MBL
Swiss	P11226 (https://www.uniprot.org/uniprotkb/P11226/entry)
别名	MBL2,COLEC1,HSMBPC,MBL,MBL2D,MBP,MBP-C,MBP1,MBPD,MBL2 Antibody,Collectin-1,Mannan-binding p rotein,Mannose-binding lectin, Mannan Binding Lectin

产品验证



Human liver was stained with anti-MBL (PT1794R) Rabbit antibody,1:200



Mouse liver was stained with anti-MBL (PT1794R) Rabbit antibody,1:200

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)